

NanKai Biotech laaja valikoima pikatestejä elintarviketuotteiden jäämien ja lisäaineiden valvontaan kattuen liha- ja siipikarja-, maito- ja meijeri-, kala- ja meri-, muna-, hunaja-, vilja- ja öljytuotteet

NKBIO Rapid Test Kits Availability:

Product Item	Product Name	Target Sample	Sensitivity (ppb)	Total Assay Time (including preparations)	Standard packing
PEN-D204R1	Penicillin Residue Rapid Test Device	Raw Milk	4 ppb	10-15 min	40/kit
β-LAC-D204R1	β-Lactam Residue Rapid Test Device	Raw Milk	3 ppb	5-10 min	40/kit
PENASE-D204R1	β-lactamase Residue Rapid Test Device	Raw Milk	3U	5-10 min	20/kit
TCs-D204R1	Tetracyclines Residue Rapid Test Device	Raw Milk	10 ppb	5-10 min	40/kit
TCs-D204R3	Tetracyclines Residue Rapid Test Device	Raw Milk	100 ppb	5-10 min	40/kit
MEL-D204R1	Melamine Residue Rapid Test Device	Raw Milk	0.2 ppm	5-10 min	40/kit
AFLM1-D204R1	Aflatoxin Residue Rapid Test Device	Raw Milk	0.1 ppb	20 min	40/kit
STR-D204R1	Streptomycin Residue Rapid Test Device	Raw Milk	50 ppb	5 min	40/kit
GEN-D204R1	Gentamicin Residue Rapid Test Device	Raw Milk	50 ppb	5-10 min	40/kit
CAP-D204R3	Chloramphenicol Residue Rapid Test Device	Raw Milk	0.3 ppb	30 min	40/kit
CAP-D204N3	Chloramphenicol Residue Rapid Test Device	Raw milk, milk powder finished Milk	0.3 ppb	30 min	40/kit
FQNs-D204R1	Fluoroquinolones Residue Rapid Test Device	Raw Milk	30 ppb	5-10 min	40/kit
SAs-D204R1	Sulfonamides Residue Rapid Test Device	Raw Milk	4 ppb	25 min	40/kit
NEO-D204R1	Neomycin Residue Rapid Test Device	Raw Milk	500 ppb	20 min	40/kit
ERY-D204R1	Erythromycin Residue Rapid Test Device	Raw Milk	40 ppb	20 min	40/kit
LIN-D204R1	Lincomycin Residue Rapid Test Device	Raw Milk	150 ppb	20 min	40/kit
TYL-D204R1	Tylosin Residue Rapid Test Device	Raw Milk	50 ppb	20 min	40/kit
KAN-D204R1	Kanamycin Residue Rapid Test Device	Raw Milk	10 ppb	10-15min	40/kit

NKBIO Rapid Test Kits Availability:

Product Item	Product Name	Sample	Sensitivity (ppb)	Total Assay Time (including preparations)	Standard packing
SAs-D203P1	Sulfonamides Residue Rapid Test Device	Tissue	4 ppb	25 min	40 /kit
SAs-D203P3	Sulfonamides Residue Rapid Test Device	Tissue	50 ppb	10 min	40/kit
FQNs-D203P3	Fluoroquinolones Residue Rapid Test Device	Tissue	10 ppb	30 min	40/kit
STR-D203P1	Streptomycin Residue Rapid Test Device	Tissue	100 ppb	10-20 min	40/kit
STR-D203P3	Streptomycin Residue Rapid Test Device	Tissue	10 ppb	30 min	40/kit
GEN-D203P1	Gentamicin Residue Rapid Test Device	Tissue	50 ppb	15-20 min	40/kit
SEM-D203P2	Furacilin Metabolite Residue Rapid Test Device	Tissue	1 ppb	120 min	40/kit
AHD-D203P2	Furadantin Metabolite Residue Rapid Test Device	Tissue	1 ppb	120 min	40/kit
AOZ-D203P2	Furazolidone Metabolite Residue Rapid Test Device	Tissue	0.5 ppb	120 min	40/kit
AMOZ-D203P2	Furaltadone Metabolite Residue Rapid Test Device	Tissue	1 ppb	120 min	40/kit
KAN-D203P1	Kanamycin Residue Rapid Test Device	Tissue	50 ppb	20 min	40/kit
ERY-D203P1	Erythromycin Residue Rapid Test Device	Tissue	50 ppb	30 min	40/kit
TYL-D203P1	Tylosin Residue Rapid Test Device	Tissue	50 ppb	30 min	40/kit
NEO-D203P1	Neomycin Residue Rapid Test Device	Tissue	500 ppb	30 min	40/kit
LIN-D203P1	Lincomycin Residue Rapid Test Device	Tissue	100 ppb	30 min	40/kit
QLAQ-D203P1	Olaquinox Residue Rapid Test Device	Tissue	4 ppb	30 min	40/kit

Veterinary residue rapid test for livestock and poultry

MRL from some main countries and areas:

* MRL data from China: Ministry of Agriculture No.235 * MRL data from EU: EU NO.508/1999/EC * MRL data from USA: USA FDA

* MRL data from Japan: JP Positive List System

Product Name	CHN MRL	EU MRL	US MRL	JP MRL
Clenbuterol	Forbidden	0.1 ppb	Forbidden	Forbidden
Ractopamine	Forbidden		50 ppb	10 ppb
Salbutamol	Forbidden			
Chloramphenicol	Forbidden	Forbidden	Forbidden	Forbidden
Nitrofurans	Forbidden	Forbidden	Forbidden	Forbidden
Fluoroquinolones	CIP+ENR 100 ppb	CIP+ENR 100 ppb	ENR 300 ppb	10-100 ppb
Tetracyclines	OTC / CTC /TC (single/total) 100 ppb	OTC / CTC /TC (single/total) 100 ppb	2 ppm	OTC / CTC /TC (total) 20 ppb DOXY 50 ppb
Sulfonamides	100 ppb	100 ppb	100 ppb	20-200 ppb
Streptomycin	600 ppb	500 ppb	500 ppb	600 ppb
Gentamycin	100 ppb	50 ppb	100 ppb	100 ppb
Erythromycin	200 ppb	200 ppb	100 ppb	50 ppb
Lincomycin	100 ppb	100 ppb	100 ppb	200 ppb
Neomycin	500 ppb	500 ppb	720 ppb	500 ppb
Olaquinox	4 ppb			
Tylosin	200 ppb	100 ppb	200 ppb	50 ppb
Chanamyn		100 ppb		100 ppb

NKBIO Rapid Test Kits Availability:

Product Item	Product Name	Target Sample	Sensitivity (ppb)	Total Assay Time (including preparations)	Standard packing
CLEN-X201U1	Clenbuterol Residue Rapid Test Device	Urine	3 ppb	3-5 min	40 /kit
CLEN-X201U2	Clenbuterol Residue Rapid Test Device	Urine	5 ppb	3-5 min	40/kits
CLEN-D201U3	Clenbuterol Residue Rapid Test Device	Urine	1 ppb	3-5 min	40/kits
RAC-X201U2	Ractopamine Residue Rapid Test Device	Urine	5 ppb	3-5 min	40/kits
RAC-D201U3	Ractopamine Residue Rapid Test Device	Urine	3 ppb	3-5 min	40/kits
SAL-X201U1	Salbutamol Residue Rapid Test Device	Urine	5 ppb	3-5 min	40/kits
SAL-D201U1	Salbutamol Residue Rapid Test Device	Urine	3 ppb	3-5 min	40/kits
CLEN-D203P2	Clenbuterol Residue Rapid Test Device	Tissue	1 ppb	40 min	40/kits
CLEN-D203P3	Clenbuterol Residue Rapid Test Device	Tissue	5 ppb	15 min	40/kits
CLEN-D203P1	Clenbuterol Residue Rapid Test Device	Tissue	3 ppb	15 min	40/kits
CLEN-X203P1	Clenbuterol Residue Rapid Test Device	Tissue	5 ppb	15 min	40/kits
RAC-D203P2	Ractopamine Residue Rapid Test Device	Tissue	1 ppb	30 min	40/kits
RAC-D203P6	Ractopamine Residue Rapid Test Device	Tissue	5 ppb	25 min	40/kits
RAC-D203P4	Ractopamine Residue Rapid Test Device	Tissue	10 ppb	15 min	40/kits
SAL-D203P1	Salbutamol Residue Rapid Test Device	Tissue	2 ppb	40 min	40/kits
SAL-D203P2	Salbutamol Residue Rapid Test Device	Tissue	5 ppb	15 min	40/kits
CAP-D203P1	Chloramphenicol Residue Rapid Test Device	Tissue	0.1 ppb	30 min	40/kits
CAP-D203P3	Chloramphenicol Residue Rapid Test Device	Tissue	0.3 ppb	30 min	40/kits
TCs-D203P1	Tetracyclines Residue Rapid Test Device	Tissue	10 ppb	20 min	40/kits
TCs-D203P2	Tetracyclines Residue Rapid Test Device	Tissue	100 ppb	20 min	40/kits